


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tryptamine hydrochloride - 99%

Product Number: 246557
CAS Number: 343-94-2
MDL: MFCD00012682
Formula: C₁₀H₁₂N₂ · HCl
Formula Weight: 196.68 g/mol



TEST

Specification

Appearance (Color)	White to Beige and Faint Brown to Brown
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃	98.5 - 101.5 %
Purity (HPLC)	≥ 98.5 %
Solubility (Turbidity)	Clear to Slightly Hazy
1% in H ₂ O	
Solubility (Color)	Colorless to Yellow and Yellow-Green

Specification: PRD.0.ZQ5.10000076548

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.