

3050 Spruce Street, Saint Louis, MO 63103, USA

 $ZrO_{2} / Y_{2}O_{3}$

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Zirconium(IV) oxide-yttria stabilized - nanopowder

 Product Number:
 572349

 CAS Number:
 308076-80-4

 MDL:
 MFCD00145546

Formula: O5Y2Zr Formula Weight: 349.03 g/mol

TEST Specification

Appearance (Color) White

Appearance (Form) Powder ICP Major Analysis Confirmed

Confirms Zr and Y Components

Surface Area (m2/g) \geq 100 Vendor Information Confirmed

Method of Manufacture Guaranteed 92% ZrO2 and 8% Y2O3 per Supplier

Specification: PRD.0.ZQ5.10000022655