

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
trans,trans-2,4-Hexadienyl acetate 97%

Product Number: W413201
CAS Number: 1516-17-2
MDL: MFCD00036521
Formula: C₈H₁₂O₂
Formula Weight: 140.18 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.472- 1.476
Infrared spectrum	Conforms to Structure
Purity (GC)	≥96.5 %
Arsenic (As)	≤3.0 ppm
Cadmium (Cd)	≤1.0 ppm
Mercury (Hg)	≤1.0 ppm
Lead (Pb)	≤10.0 ppm

Remarks:

Specification Date : 06/22/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.