


Product Specification

Product Name:
(±)-Miconazole nitrate salt - imidazole antibiotic

Product Number: **M3512**
 CAS Number: 22832-87-7
 Formula: C₁₈H₁₄Cl₄N₂O·HNO₃
 Formula Weight: 479.1 AMU



TEST	Specification
Appearance (Form)	Powder
Appearance (Color)	White to Off White
Specific Rotation (BP)	-0.10 - 0.10 °
Purity (HPLC) > = 98.0%	Conforms
Identification IR spectrum conforms to control or reference standard	Conforms
Expiration Date Period 4 Years	-----

Specification: PRE.1.ZQ5.10000010100

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

