

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Monoclonal Anti-c-Myc antibody produced in mouse – clone 9E10, purified immunoglobulin, buffered aqueous solution

Product Number: **M4439**

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	2.0 - 3.0 mg/ml
WB-Cell Line/Tissue Extract	293-T Cells or COS7 or E.coli Expressing c-Myc Tag
WB Titer	1:5,000-1:7,000
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000015073

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.