

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-[[4'-(Trifluoromethoxy)phenoxy]methyl]phenylboronic acid

Product Number:

652156

CAS Number:

849062-07-3

MDL:


MFCD07369751

Formula:

C₁₄H₁₂BF₃O₄

Formula Weight:

312.05 g/mol



TEST

Specification

Appearance (Color)

White to Off White

Appearance (Form)

Conforms to Requirements

Powder, Crystals or Crystalline Powder, and/or Chunks

Proton NMR Spectrum

Conforms to Structure

Titration with NaOH

≥ 95.0 %

Specification: PRD.0.ZQ5.10000039763

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.