


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
sSPhos

Product Number: 677280
CAS Number: 1049726-96-6
MDL: MFCD07781994
Formula: C₂₆H₃₄NaO₅PS · xH₂O
Formula Weight: 512.57 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Yellow to Brown
Appearance (Form)	Conforms to Requirements
Solid, Powder or Crystals	
Proton NMR Spectrum	Conforms to Structure
Phosphorus NMR Spectrum	Conforms to Structure
Carbon (dry basis)	58.7 - 63.1 %
Water (by Karl Fischer)	1.7 - 5.0 %

Specification: PRD.3.ZQ5.10000018774

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

