

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Phenol:Chloroform 5:1 - for molecular biology

Product Number:

P1944

Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form) Mixture	Liquid
% Purity	≥ 99.0
UV Absorbance 330nm	≤ 0.2
UV Absorbance 405nm	≤ 0.1
UV Absorbance 510nm	≤ 0.1
pH	4.3 - 4.7

Specification: PRD.0.ZQ5.10000049019

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.