


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,5-Cyclooctadiene - contains 50-150 ppm 4-tert-Butylcatechol, 99%

Product Number: C109207
CAS Number: 111-78-4
MDL: MFCD00001752
Formula: C₈H₁₂
Formula Weight: 108.18 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %
Stabilizer	50 - 150 ppm
4-tert-Butylcatechol	

Specification: PRD.1.ZQ5.10000019374

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.