

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Nickel(II) chloride hexahydrate - ReagentPlus®

Product Number:**223387**NiCl₂ • 6H₂O

CAS Number:

7791-20-0

MDL:

MFCD00149809

Formula:

Cl₂Ni • 6H₂O

Formula Weight:

237.69 g/mol

TEST**Specification**

Appearance (Color)

Conforms to Requirements

Green

Appearance (Form)

Powder or Crystals

ICP Major Analysis

Confirmed

Confirms Nickel Component

Compleximetric EDTA

23.0 - 26.0 %

% Nickel

Recommended Retest Period

36 Months

Specification: PRD.5.ZQ5.10000005655

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.