


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium thiosulfate - free-flowing, Redi-Dri™, ReagentPlus®, 99%

Product Number: 909173
CAS Number: 7772-98-7



Formula: $\text{Na}_2\text{S}_2\text{O}_3$
Formula Weight: 158.10 g/mol

TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder or Crystals or Granules
Titration by Iodine	98.5 - 101.5 %
ICP Major Analysis	Confirmed
Confirms Sodium and Sulfur Components	
Solubility (Turbidity)	Clear
100 mg/ml (10%), H ₂ O	
Solubility (Color)	Colorless
Recommended Retest Period	-----
3 Years	

Specification: PRD.2.ZQ5.10000098895

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

