

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tricyclo[5.2.1.0^{2,6}]decanedimethanol diacrylate

Product Number: 496669
CAS Number: 42594-17-2



Formula: C₁₈H₂₄O₄
Formula Weight: 304.38 g/mol
Storage Temperature: 2 - 8 °C

| TEST | Specification |
|--------------------------------------|-----------------------|
| Appearance (Color) | Colorless |
| Appearance (Form) | Liquid |
| Refractive index at 20 °C | 1.503 - 1.509 |
| Proton NMR spectra | Conforms to Structure |
| Specification: PRD.1.ZQ5.10000017646 | |

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.