


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,6-Difluorostyrene - 99%, contains 0.25% 4-tert-butylcatechol as inhibitor

Product Number: **374407**
CAS Number: 207226-37-7
MDL: MFCD00134389
Formula: C₈H₆F₂
Formula Weight: 140.13 g/mol



TEST

Specification

Appearance (Color)	Colorless to Yellow
Sediment may form on standing, but does not affect use	
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %

Specification: PRD.0.ZQ5.10000014256

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.