


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
trans-2-Hexen-1-al - natural, ≥95%, FG

Product Number: **W256009**
CAS Number: 6728-26-3



Formula Weight: 98.15 g/mol
Storage Temperature: 2 - 8 °C

TEST

Specification

Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Refractive index at 20 °C	1.445 - 1.449
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 95.00 %
Organoleptics	Pass
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
60 Months	

Specification: PRD.1.ZQ5.1000060396

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.