

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification


Product Name:

Chitosan - medium molecular weight

 Product Number:
 448877

 CAS Number:
 9012-76-4

 MDL:
 MFCD00161512



TEST

Appearance (Color)

Appearance (Form)

Powder and/or Chips

Deacetylation

Viscosity

c = 1%, 1% Acetic Acid

Specification: PRD.1.ZQ5.10000030084

Specification

Off White to Beige and Faint Brown to Light Brown

Conforms to Requirements

75 - 85 %

200 - 800 cps

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.