

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


11-Mercaptoundecylphosphonic acid - 95%

Product Number:

754269

CAS Number:

156125-36-9



Formula:

C₁₁H₂₅O₃PS

Formula Weight:

268.35 g/mol

Storage Temperature:

2 - 8 °C

TEST

Specification

Appearance (Color)

White to Yellow to Beige

Appearance (Form)

Powder or Crystals

Proton NMR spectrum

Conforms to Structure

Phosphorus NMR Spectra

Conforms to Structure

Carbon

46.5 - 51.9 %

Sulfur

11.3 - 12.6 %

Specification: PRD.1.ZQ5.10000036741

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.