

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Pectolyase from *Aspergillus japonicus* – plant cell culture tested, lyophilized powderProduct Number: **P5936**

MDL: MFCD00131813

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Activity	≥ 0.3
Units/mg solid	
One unit will liberate 1.0 micromole galacturonic acid from polygalacturonic acid per min, pH 5.5, 25 deg C	
Plant Cell Culture Test	Pass

Specification: PRD.4.ZQ5.10000007797

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.