

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Copper(I) bromide - 99.999% trace metals basis

Product Number: 254185
CAS Number: 7787-70-4
MDL: MFCD00010969
Formula: BrCu
Formula Weight: 143.45 g/mol

CuBr

TEST	Specification
------	---------------

Appearance (Color)	Conforms to Requirements
Pale Green	
Appearance (Form)	Powder
Titration with EDTA	43.6 - 45.0 %
% Cu	
ICP Major Analysis	Confirmed
Confirms Copper Component	
Purity	Conforms
99.999% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 15.0 ppm

Specification: PRD.1.ZQ5.10000046190

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.