


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
N-[3-(Dimethylamino)propyl]methacrylamide - 99%, contains MEHQ as inhibitor

Product Number: 409472
CAS Number: 5205-93-6
MDL: MFCD00038359
Formula: C₉H₁₈N₂O
Formula Weight: 170.25 g/mol



TEST	Specification
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Viscous Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %
MEHQ as Inhibitor	600 - 750 ppm

Specification: PRD.1.ZQ5.10000005789

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.