

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ethyl (S)-(-)-4-chloro-3-hydroxybutyrate - 96%

Product Number:

460524

CAS Number:

86728-85-0

MDL:


MFCD00211241

Formula:

C₆H₁₁ClO₃

Formula Weight:

166.6 g/mol



TEST

Specification

Appearance (Color)

Colorless to Light Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 95.5 %

Enantiomeric Excess

≥ 96.5 % EE

Specification: PRD.0.ZQ5.10000037698

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.