

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(Iodomethyl)triphenylphosphonium iodide - 95%

Product Number: **725285**
CAS Number: 3020-28-8

Formula: C₁₉H₁₇I₂P
Formula Weight: 530.12 g/mol



TEST

Specification

Appearance (Color)	White to Yellow to Beige
Appearance (Form)	Powder or Solid or Crystals
Proton NMR spectrum	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure
Titration with AgNO ₃ After Oxygen Combustion	94.5 - 105.5 %

Specification: PRD.1.ZQ5.10000050928

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.