

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Bis(3,5-di(trifluoromethyl)phenyl)chlorophosphine

Product Number:

694746

CAS Number:

142421-57-6

MDL:


MFCD01630852

Formula:

C₁₆H₆ClF₁₂P

Formula Weight:

492.63 g/mol

**TEST****Specification**

Appearance (Color)

Colorless and White to Off White

Appearance (Form)

Liquid or Solid

Proton NMR spectrum

Conforms to Structure

Phosphorus NMR Spectra

Conforms to Structure

Specification: PRD.0.ZQ5.10000052843

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.