


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1-(2,3-Dichlorophenyl)piperazine hydrochloride - 97%

Product Number: 679135
CAS Number: 119532-26-2
MDL: MFCD00190238
Formula: C₁₀H₁₃Cl₂N₂
Formula Weight: 267.58 g/mol



TEST

Specification

Appearance (Color)	White to Tan
Appearance (Form)	Crystal to Powder
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃	96.5 - 103.5 %
Purity (HPLC)	≥ 96.5 %

Specification: PRD.1.ZQ5.10000016189

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.