

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Elafin human – >90% (by MS, HPLC and SDS-PAGE), recombinant, expressed in *Saccharomyces cerevisiae*Product Number: **E7280**

MDL: MFCD00146824

Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	White to Faint Yellow
Appearance (Form)	Powder
Purity	> 90 %
by SDS Electrophoresis	
Purity (HPLC)	> 90 %
Activity	> 90 %
Inhibition of Human Leukocyte Elastase (%)	

Specification: PRD.1.ZQ5.10000007111

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.