

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-Ethylhexyl trans-4-methoxycinnamate - 98%, contains 500-1000 ppm BHT as stabilizer

Product Number:

437174

CAS Number:

83834-59-7

MDL:


MFCD00072582

Formula:

C₁₈H₂₆O₃

Formula Weight:

290.40 g/mol

**TEST****Specification**

Appearance (Color)

Colorless to Light Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 97.5 %

BHT as Stabilizer

Confirmed

500-1000 ppm

Specification: PRD.1.ZQ5.10000031319

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.