

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

B-Iodo-9-BBN solution - 1.0 M in hexanes

Product Number:

416908

MDL:


MFCD00192226

Formula:

C₈H₁₄BI

Formula Weight:

247,91 g/mol



TEST

Specification

Appearance (Color)

Colorless to Very Dark Yellow to Very

Appearance (Form)

Dark Brown and Colorless to Very Dark Red to Very Dark Purple

Titration with AgNO₃

Liquid

0.95 - 1.10 M

Specification: PRD.3.ZQ5.10000047186

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.