

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methyl acetimidate hydrochloride - technical grade

Product Number: 254940
CAS Number: 14777-27-6
MDL: MFCD00012571
Formula: C₃H₇NO · HCl
Formula Weight: 109.55 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form) Powder, Crystals, Crystalline Powder and/or Chunks	Conforms
Proton NMR spectra MAY CONTAIN UP TO 15% NH ₄ CL	Conforms to Structure
Titration by AgNO ₃	88.0 - 112.0 %

Specification: PRD.0.ZQ5.10000052882

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.