


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tris(2-chloroethyl) phosphate - 97%

Product Number: 119660
CAS Number: 115-96-8
MDL: MFCD00000967
Formula: C₆H₁₂Cl₃O₄P
Formula Weight: 285.49 g/mol



TEST	Specification
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Viscous Liquid
Infrared spectrum	Conforms to Structure
Titration with AgNO ₃ After Oxygen Combustion	96.5 - 103.5 %

Specification: PRD.1.ZQ5.10000040957

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.