

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Butyl sulfide - ≥95%, FG

Product Number: W221503
CAS Number: 544-40-1
MDL: MFCD00009468
Formula: C8H18S
Formula Weight: 146.29 g/mol



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.450 - 1.454
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 95.0 %
Specific Gravity (at 25 deg C)	0.820 - 0.860
Color test (APHA) 200 Max	Pass
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period 5 Years	-----

Specification: PRD.4.ZQ5.10000037109

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.