

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-phospho-DAP-Kinase (pSer³⁰⁸) antibody produced in mouse – clone DKPS308, purified immunoglobulin, buffered aqueous solutionProduct Number: **D4941**

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	2.0 - 2.5 mg/ml
WB-Cell Line/Tissue Extract	293-T Cells Exp. DAP Kinase
WB Working Concentration	1 - 2 µg/ml
Recommended Retest Period	-----
4 Years	

Specification: PRE.1.ZQ5.10000016485

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.