


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2'-Hydroxyflavanone - ≥98%

Product Number: H3780
CAS Number: 17348-76-4
MDL: MFCD00017702
Formula: C₁₅H₁₂O₃
Formula Weight: 240.25 g/mol



TEST	Specification
Appearance (Color)	White to Faint Yellow
Appearance (Form)	Powder
Solubility (Color)	Faint Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
50 mg/mL MeOH	
Proton NMR spectrum	Conforms
¹³ C NMR Spectrum	Conforms
Purity (TLC)	≥ 98 %

Specification: PRD.1.ZQ5.10000074874

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.