

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Yeast from *Saccharomyces cerevisiae* - Type II

Product Number:

YSC2

Storage Temperature:

2 - 8 °C

### TEST

### Specification

Appearance (Color)

Faint Beige to Brown

Appearance (Form)

Powder or Solid

Activity

Conforms

Enzymatic Activity: Contains 3-Phosphoglycerate Kinase Activity.

Unit Definition: One unit will convert 1.0 micromole

of 1,3-Diphosphoglycerate to 3-Phosphoglycerate per minute

at pH 6.9 at 25 deg C

Specification: PRD.0.ZQ5.10000043486

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.