

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)

Email USA: [techserv@sial.com](mailto:techserv@sial.com)


Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
FK 269 Co(III) PF<sub>6</sub> salt

Product Number: **805521**

Formula: C<sub>22</sub>H<sub>18</sub>CoF<sub>18</sub>N<sub>10</sub>P<sub>3</sub>  
Formula Weight: 916.27 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Faint Orange to Very Dark Orange
--------------------	----------------------------------

Appearance (Form)	Powder or Crystals
-------------------	--------------------

Infrared Spectrum	Conforms to Structure
-------------------	-----------------------

Carbon Content	29.4 - 31.0 %
----------------	---------------

Nitrogen Content	15.6 - 16.4 %
------------------	---------------

Product of Supplier	Confirmed
---------------------	-----------

Greatcell Solar®

Registered Trademark

Greatcell Solar is a registered trademark of Greatcell Solar

Recommended Retest Period

2 Years

Note

-----

Specification: PRD.0.ZQ5.10000099214

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.