

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Carbon nanotube, multi-walled - carboxylic acid functionalized, thin, extent of labeling: >8% carboxylic acid functionalized, avg. diam. \times L 9.5 nm \times 1.5 μ m

Product Number: 755125

Formula: C



TEST	Specification
Appearance (Color)	Black
Appearance (Form)	Powder
Carbon Content (TGA)	> 80 %
Miscellaneous Assay -COOH Functionalization (XPS)	> 8
Product of Supplier Nanocyl Inc	Confirmed

Specification: PRD.1.ZQ5.10000028402

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.