

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Activated charcoal – 4–14 mesh, granular, Norit® PK 3–5

Product Number:**292591****C**

CAS Number:

7440-44-0

MDL:

MFCD00133992

Formula:

C

Formula Weight:

12.01 g/mol

TEST**Specification**

Appearance (Color)

Black

Appearance (Form)

Granules

Size

Confirmed

4 - 14 Mesh

Registered Trademark

Confirmed

Norit is a registered trademark of Norit N.V.

Specification: PRD.0.ZQ5.10000022476

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.