


Product Specification

Product Name:
Apicidin - ≥98% (HPLC), from microbial

Product Number: **A8851**

Formula: C₃₄H₄₉N₅O₆

Formula Weight: 623.78 AMU



TEST	Specification
Appearance (Form)	Powder
Appearance (Colour)	White to White w/ Yellow Cast
Solubility (Solvent)	DMSO
Solubility (Conc)	0.90 - 1.10 mg/ml
Solubility (Heating)	Yes
Solubility (Sonication)	Yes or No
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless
Purity (HPLC)	≥ 98.00 %
NMR (Solvent)	CDCl ₃
Identity by NMR	Consistent
Recommended Retest Period	-----
2 Years	

Specification: PRE.3.ZQ5.10000002005

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

