

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

trans,trans-2,4-Hexadienoic acid potassium salt - ReagentPlus®, 99%

Product Number:

359769

CAS Number:

24634-61-5

MDL:


MFCD00016546

Formula:

C₆H₇KO₂

Formula Weight:

150.22 g/mol



TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Conforms to Requirements

Granules or Beads

Infrared spectrum

Conforms to Structure

Purity (Titration by HCL04)

98.5 - 101.0 %

Specification: PRD.0.ZQ5.10000053790

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.