

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Desmin antibody produced in mouse – clone DE-U-10, ascites fluid

Product Number: D1033

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Very Faint Yellow to Very Dark Yellow to Very Dark Brown
Protein by Biuret	15 - 55 mg/ml
Mouse IgG1 by RID	6.0 - 10.0 mg/ml
WB-Cell Line/Tissue Extract	Chicken Gizzard/C2C12
WB Titer	1:200-1:400
IHC - Tissue	Human Appendix
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:40-1:80
Recommended Retest Period	-----
5 Years	

Specification: PRE.8.ZQ5.10000015671

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

