

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-O- α -D-Glucopyranosyl-L-ascorbic Acid - $\geq 98\%$ (HPLC)

Product Number:

SMB00390

CAS Number:

129499-78-1

Formula:


C₁₂H₁₈O₁₁

Formula Weight:

338.26 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

White to Off White

Appearance (Form)

Powder

Proton NMR Spectrum

Conforms to Structure

Purity (HPLC)

 $\geq 98\%$

Specification: PRD.0.ZQ5.10000067176

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.