

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)

Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

**Product Specification**

Product Name:

(2R,2'R-(+)-[N,N'-Bis(2-pyridylmethyl)]-2,2'-bipyrrolidine)bis(acetonitrile)iron(II) hexafluoroantimonate

Product Number:


**731382**

Formula:

C<sub>24</sub>H<sub>32</sub>F<sub>12</sub>FeN<sub>6</sub>Sb<sub>2</sub>

Formula Weight:

931.90 g/mol



**TEST**

**Specification**

Appearance (Color)	Conforms to Requirements
Off White to Brown and Off White to Purple	
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	29.8- 32.0 %
Nitrogen	7.9- 10.1 %
Optical Rotation	20.0- 30.0 deg
c = 0.2%; Acetonitrile	
ICP Major Analysis	Confirmed
Confirms Iron and Antimony Components	

**Remarks:**

Specification Date : 11/30/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.