

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium hydride - 60 % dispersion in mineral oil

Product Number: 452912

MDL: MFCD00003471
Formula: HNa
Formula Weight: 24.00 g/mol

NaH

TEST	Specification
------	---------------

Appearance (Color)	Conforms to Requirements
Faint Violet to Grey to Blue-Grey	
Appearance (Form)	Conforms to Requirements
Moist, Paste-Like Powder or Dispersion	
Hydrogen Evolution	57.0 - 63.0 %
Dispersion	
Expiration Date Period	-----
3 Years	

Specification: PRD.2.ZQ5.10000004904

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.