


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Hydrazino-2-imidazoline hydrobromide - 98%

Product Number: 197173
CAS Number: 55959-84-7
MDL: MFCD00012694
Formula: C₃H₈N₄ · HBr
Formula Weight: 181.03 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Assay	≥ 97.5 %

Specification: PRD.1.ZQ5.10000022728

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

