

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,3-Bis(trifluoromethyl)-5-bromobenzene - 99%

Product Number:

290157

CAS Number:

328-70-1

MDL:


MFCD00000381

Formula:

C₈H₃BrF₆

Formula Weight:

293.00 g/mol

**TEST****Specification**

Appearance (Color)

Conforms to Requirements

Colorless to Light Pink to Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 99 %

Specification: PRD.0.ZQ5.10000037637

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.