

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hyaluronic acid sodium salt from *Streptococcus zooepidemicus* - bacterial glycosaminoglycan polysaccharide

Product Number:


H9390

CAS Number:

9067-32-7

MDL:

MFCD01773053



Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Color)

White

Appearance (Form)

Powder

Solubility (Color)

Colorless

Solubility (Turbidity)

Clear

at 10 mg plus 0.5 mL of H₂O

Suitability

Suitable

Suitable for use as substrate for hyaluronidase

Recommended Retest Period

3 years

Specification: PRD.0.ZQ5.10000047458

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.