


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium bis(trimethylsilyl)amide - 95%

Product Number: 235083
CAS Number: 1070-89-9
MDL: MFCD00009835
Formula: C₆H₁₈NNaSi₂
Formula Weight: 183.37 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Faint Yellow to Faint Beige
Appearance (Form)	Powder or Chunks
Infrared spectrum	Conforms to Structure
Titration by HCL	94.0 - 106.0 %

Specification: PRD.1.ZQ5.10000026216

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.