

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Aluminum – pellets, 3–8 mesh, ≥99.999% trace metals basis

Product Number: 326941
CAS Number: 7429-90-5
MDL: MFCD00134029
Formula: Al
Formula Weight: 26.98 g/mol

Al

TEST	Specification
Appearance (Color)	Conforms to Requirements
Silver	
Appearance (Form)	Pellets
Size	Confirmed
3-8 Mesh	
ICP Major Analysis	Confirmed
Confirms Aluminum Component	
Trace Metal Analysis	≤ 10.0 ppm
Purity	Meets Requirements
≥99.999% Based on Trace Metals Analysis	

Specification: PRD.0.ZQ5.10000017538

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.