

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Benzoin – purified by sublimation, ≥99.5%

Product Number:

399396

CAS Number:

119-53-9

MDL:


MFCD00004496

Formula:

C₁₄H₁₂O₂

Formula Weight:

212.24 g/mol

**TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Carbon

78.8 - 79.6 %

Purity (HPLC)

≥ 99.5 %

Specification: PRD.0.ZQ5.10000026367

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.