


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Podocarpic acid - 98%

Product Number: 119792
CAS Number: 5947-49-9
MDL: MFCD00074825
Formula: C₁₇H₂₂O₃
Formula Weight: 274.35 g/mol



TEST

Specification

Appearance (Color)	White to Tan
Appearance (Form)	Conforms to Requirements
Powder, Crystals or Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Titration with NaOH	97.5 - 102.5 %
Optical Rotation	129.0 - 137.0 deg
C = 4%, 95% Ethanol	

Specification: PRD.2.ZQ5.10000030141

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.