

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Lithium hexafluorophosphate - 98%

**Product Number:** 201146  
**CAS Number:** 21324-40-3  
**MDL:** MFCD00011096  
**Formula:** F<sub>6</sub>LiP  
**Formula Weight:** 151.91 g/mol

LiPF<sub>6</sub>

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Chunks
ICP Major Analysis Confirms Lithium and Phosphorus Components	Confirmed
Gravimetric Analysis With Tetraphenylarsonium Chloride	97.5 - 102.5 %

Specification: PRD.1.ZQ5.10000032477

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.