

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium phosphate tribasic – reagent grade, ≥98%

Product Number:

P5629

CAS Number:

7778-53-2

MDL:

MFCD00036295

Formula:

K₃O₄P

Formula Weight:

212,27 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/ml, H ₂ O	
Loss on Ignition	≤ 5.0 %
Titration with HCl	≥ 98.0 %
Ignited basis	
Iron (Fe)	≤ 0.001 %
Recommended Retest Period	-----
4 Years	

Specification: PRD.1.ZQ5.10000016298

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.